



# DNA Extraction, Gel Electrophoresis & Observation of Polytene Chromosomes



# DNA Fingerprinting

Deoxyribonucleic acid (DNA) is the genetic material that acts as the blueprint for protein synthesis by cells. In mammals, a large fraction of the total DNA does not code for proteins. Polymorphic DNA refers to chromosomal regions that vary among individuals. By examining several of these regions within genomic DNA, one can determine a “DNA Fingerprint” for an individual. DNA polymorphisms are now widely used for determining paternity/ maternity, identification of human remains and to determine the genetic basis of various inherited diseases. The most widely used and far-reaching application has been to field of criminal forensics. DNA from crime victims and offenders can now be definitively matched, affecting outcomes of criminal and civil trials.

## Equipment and Materials

250 ml conical flask .....	1
Cooling basin.....	1
Digital thermometer .....	1
Electronic balance .....	Share
Gel bed with comb (Close off the open ends by paper tape) .....	1 set
Microwave oven.....	Share
Microwave oven wrap.....	Share
Latex gloves .....	1 pair
Pipettement and tips.....	1 set

## Chemicals

Agarose .....	Share
TAE buffer(1X).....	50ml

# Gel Electrophoresis

## Procedures

### Preparation of Agarose Gel (1.5% , 7 cm x 10 cm Gel)

1. Put the conical flask on the electronic balance and set Zero.
2. Add 0.75g agarose into the conical flask.
3. Add 50 ml TAE buffer into the conical flask.
4. Wrap the opening of conical flask and make a hole on it.
5. Use a microwave oven to heat the mixture until all the agarose is completely dissolved. (Clear solution can be observed)
6. Cool the agarose solution to 55 °C with a cooling basin.
7. Add 5µl SYBR red stain. (Done by teacher)
7. Pour the cooled agarose solution into the gel bed.
8. Allow the gel to completely solidify. It will become firm after approximately 20 minutes.
9. Remove the comb by slowly pulling straight up. Be careful and evenly to prevent tearing the sample wells.



### Load the samples and Run the Gel

1. Load 25µl of the DNA samples from tubes A-E into the wells. The content of Tubes A-E is shown as follow,
  - A. DNA Standard marker
  - B. Child's DNA fingerprint
  - C. Mother's DNA fingerprint
  - D. Male 1 DNA fingerprint
  - E. Male 2 DNA fingerprint
2. After DNA samples are loaded, connect the apparatus to the D.C. power source and set the power source at the required voltage.
3. Check that current is flowing properly – you should see bubbles forming on the two platinum electrodes. Conduct electrophoresis for the length of time specified by your instructor.
4. After electrophoresis is completed, proceed to DNA visualization by exposing the gel in UV light

## DNA extraction from strawberries

### Equipment and Materials

1ml Syringe.....	1
5ml Syringe.....	1
100ml Beaker .....	1
Coffee filter .....	1
Filter funnel.....	1
Pestle and Mortar .....	1 set
Fruit Sample	
Test tubes .....	2

### Chemicals

Extraction Buffer .....	10 ml
100 ml detergent	
900 ml distilled water	
15 grams NaCl	
Cold Alcohol(pre-measured)...	4 ml

### Procedures


#### DNA Isolation using Strawberries (wear eye protection, or be careful).

1. Mash the strawberries with pestle and mortar for 2 minutes.
2. Add 10 ml DNA Extraction buffer (Pre-measured).
3. Mash for 1 minute.
4. Filter solution by coffee filter into beaker
5. Transfer 2 ml of this solution back into the two test tubes. (It will be very red)
6. Add 4 ml (2X the volume of strawberry solution) of 75% cold ethanol and 98% cold ethanol to the test tubes respectively. (Add the ethanol slowly so that it stays on top of the strawberry solution.
7. Record the observations


### NOTE

If the tube just sits, DNA will precipitate at the interface with bubbles. More DNA can be observed if the tube was shake in gently and eventually float to the top.


## Principles




① DNA is in the nucleus of the cell



② Cell membrane is disrupted with a detergent



③ Cold Alcohol is added to the tube slowly to separate DNA from other cell components. DNA moves to the alcohol layer.



④ DNA is spooled onto wooden sticks.

## Discussion

1. What is the function of extraction buffer?
2. What is the function of alcohol? Why it should be cold?
3. Why the fruit sample should be frozen beforehand?
4. Why we usually select strawberries and Kiwi fruit as sample for DNA extraction?
5. Why the yields of DNA extraction are different within different fruit sample?
6. How to increase the yield of DNA extraction? Why?
7. Is the precipitate as a pure DNA?

# Observation of Polytene Chromosomes in Larva of *Chironomidae*

Polytene chromosomes are over-sized chromosomes, and are commonly found in cells which undergo repeated rounds of DNA replication without cell division, called endocycle. As the sister chromatids remain synapsed together, polytene chromosomes will be formed after multiple rounds of replication occurs. This kind of specialized cells may be found in the salivary glands of *Drosophila*.



## Equipment and Materials

Chironomidae Larva ..... 1 pack  
Forceps.....2  
Filter paper .....few sheets  
Slide.....3  
Coverslip.....3  
Stereomicroscope..... 1  
Light Microscope ..... 1

## Chemicals

Aceto-orcein ..... Share  
Saline..... Share

## Procedures

1. Prepare a clean slide with 2-3 drops of saline.
2. Select a large larva and place it on the slide.
3. Put the slide on a stereomicroscope.
4. Use pins or forceps to grasp the larva by its middle part and its head respectively.
5. Gently stretch the larva by pulling on it until its head separates from the rest of its body.
6. Look for the salivary glands in the head section. The glands are very small, fairly transparent, usually paired and have dark fat particles attached.
7. Soak the saline with filter paper.
8. Add 2-3 drops of aceto-orcein stain to the salivary glands and keep it for 5-10 mins.
9. After the stain has set, cover the salivary glands with coverslip.
10. Place your thumb on the over the coverslip and press down slowly.

## Discussion

1. What are the benefits of having polyteny in the salivary glands?
2. What are the similarities & differences between polyteny and polyploidy?